boso medistar S

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2

**User
Instructions**

35

**Mode
d'emploi**

68

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101

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134

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

Table of contents

Product features	37
Contents of package	38
Quick guide	39
Introduction	40
Blood pressure values	42
Systolic and diastolic blood pressure value	42
WHO blood pressure cut-off values	42
Starting up	44
Inserting batteries	44
Battery status	45
Changing batteries	45
General instructions for self-measurement	46
Preparing to measure your blood pressure	48
Note on cuff sizes	48
How to wear the unit	49

36 Table of contents


Measuring your blood pressure	52
Interrupting a measurement	53
Measurement display	54
Completing the blood pressure record card / Retrieving the average value	56
Displaying individual stored values	57
Error messages	58
Cleaning and Disinfection	60
Warranty	60
Customer service / Disposal	61
Technical data	62
Calibration checks - Testing instructions (for trained specialist staff only)	65

Product features

- | | | |
|---|--|-----------------|
|  | Battery status display | see page 45 |
|  | Measurement in progress | see page 53 |
|  | Irregular heartbeat during measurement | see page 55, 59 |
|  | Stored value | see page 57 |

38 Contents of package

1 Blood pressure monitor
boso-medistar S

1 Warranty certificate


1 Case

2 LR 03 (AAA)

1 Blood pressure record card

1 User instructions

Gebruiksaanweisung	2
User Instructions	35
Mode d'emploi	68
Manuale di istruzioni	101
Manual del Usuario	134

Quick guide

39

➔ Insert batteries, ensure correct positioning (page 44)

➔ Attach unit (page 49)

 When measuring, the wrist must be held at the same height as the heart!

➔ Start (and stop) a measurement by pressing the 'START' button (page 52)

Measured systolic, diastolic and pulse values are displayed after the measurement is completed (page 54)

➔ Retrieving values from memory (page 56): with the instrument switched off, press the 'MEMORY' but-

ton. The first value displayed is the average value of all stored measurements. The total number of stored measurements is displayed at the same time, e.g. 'A30'. After repeated pressing the MEMORY button, each individual measurement is displayed.

 This quick guide does not replace the detailed information on the use and safety of the blood pressure monitor featured in the user instructions.

Therefore, please ensure reading the user instructions!

40 Introduction

Dear Customer,

We are delighted that you have decided to purchase a boso blood pressure monitor. The boso brand is synonymous with optimal quality and precision and is also number 1 among professionals: 77% of all German doctors currently use boso blood pressure monitors (GfK 1/2010).

The experience gained from decades of use by professionals has also been incorporated into all the patient monitors for home measurement.

This instrument has also passed our strict quality control procedures and is

your reliable partner for monitoring your blood pressure values.

 Please read these user instructions before using the instrument for the first time, as correct blood pressure readings can only be obtained if the instrument is operated correctly.

In these user instructions, '➡' indicates an action by the user.

To get help concerning usage or maintenance please contact your dealer or the manufacturer (contact details can be found on the back cover of this manual).

Introduction

41


An unexpected operating condition or an incidence which has worsened the health condition or could have worsened the health condition should be reported to the manufacturer immediately.

If the instrument is to be sold, ensure that these user instructions are enclosed.

If the monitor is used in medical practice (in accordance with the legislation governing the operation of medical devices), it should be calibrated at regular intervals (see page 65).

The boso-medistar S blood pressure monitor is suitable for patients of all ages with a wrist circumference of between 13,5 and 19,5 cm. The monitor is not suitable for neonatal.

Wireless communication devices, such as home networking devices, mobile phones, cordless phones and their base stations, walkie-talkies can affect this blood pressure monitor. Therefore, a minimum distance of 11 feet should be kept from such devices.

42 Blood pressure values


To determine the blood pressure, **two values** need to be measured:

- **Systolic (upper) blood pressure:**
This is produced when the heart muscle contracts and blood is forced into the blood vessels.

- **Diastolic (lower) blood pressure:**
This is present when the heart muscle expands and fills with blood again.

Blood pressure measurement values are expressed in mmHg (mm of mercury).

The World Health Organization (WHO, 1999) has defined the following cut-off values for assessing blood pressure:		
	systolic	diastolic
High blood pressure	over 140 mmHg	over 90 mmHg
Borderline normal blood pressure	130 to 139 mmHg	85 to 89 mmHg
Normal blood pressure	120 to 129 mmHg	80 to 84 mmHg
Optimum blood pressure	up to 119 mmHg	up to 79 mmHg

Blood pressure values

43


Whether medical treatment is necessary depends not only on the blood pressure, but also on the patient's risk profile. Please contact your GP if one of the values (systole, diastole) consistently exceeds the cut-off point for high blood pressure.

Self-measurement of blood pressure often produces somewhat lower values at home than at the doctor's.

The German Hypertension League has therefore defined lower cut-off values for home measurements:

Home measurement: 135/85 mmHg

Measurement in the
doctor's surgery: 140/90 mmHg

44 Starting up

The device must be installed and put into operation in accordance with the instructions in this manual.


Illustration 1

⚠ Never mix old and new batteries or different makes.

⚠ If the batteries are not inserted correctly, the monitor will not work and this may cause the batteries to heat up and leak, damaging the monitor.

➔ Inserting the batteries

The battery compartment is located on the left side of the monitor. Insert batteries in the battery compartment as shown in illustration 1.

➔ Use only high quality leak-proof batteries of the correct specifications (see 'Technical data' page 62).

➔ Remove the batteries if the monitor is not to be used for any length of time.

Starting up

45

The monitor has a battery status display (see page 37)

No further measurements possible.
Replace battery.

⚠ Change of battery deletes memory!

The removal of batteries will erase the memory store immediately.


⚠ Care for the environment

Do not dispose of used batteries and charger units in domestic waste. You can take these to collection sites for used batteries or special waste. Contact your local authority for information.

46 General instructions for self-measurement

- 
1. Fluctuations in blood pressure are quite normal. Even when measurements are repeated, there can be marked differences. Single or irregular measurements do not provide reliable information about the actual blood pressure. A reliable evaluation is only possible if regular measurements are taken under comparable conditions and the measurement values are entered in the blood pressure record card.
 2.  Self-testing does not mean self-treatment. Do not on **any** account change the dose of medicine prescribed by your doctor of your own accord.
 3. Irregular heartbeats can affect the accuracy of the monitor or result in incorrect measurements (see page 59).
 4. Incorrect measurements can also occur in patients with a weak pulse (possibly in patients wearing a pacemaker, for example). The blood pressure monitor itself has no effect on the pacemaker.
 5. If you are pregnant, consult your doctor before using the monitor.

General instructions for self-measurement


47

6. Always take the measurement in a calm and relaxed state.

You are recommended to check your blood pressure twice a day – in the morning on getting up and in the evening once you have relaxed after work.

7. Always measure your blood pressure from the wrist with the higher blood pressure values (where not required otherwise (see page 50)).

To do this, measure the blood pressure first in both wrists and thereafter always from the wrist with the higher blood pressure.

48 Preparing to measure your blood pressure

Note on cuff sizes:

The unit must only be used in combination with the CW 42 cuff (wrist circumference 13.5 to 19.5 cm).

GB


Observe the following instructions when measuring your blood pressure:

Avoid smoking or drinking coffee for one hour before measuring your blood pressure.

Sit comfortably for the blood pressure measurement. Support your back and arms. Do not cross your legs. Place your feet flat on the floor. Relax for 5 minutes before taking a measurement. Do not move during the measurement.

Preparing to measure your blood pressure

49

How to wear the unit

The measurement should be taken from the bare wrist.

⚠ Please make sure the display is positioned at your inner wrist (see illustration 2).


Illustration 2

➔ Close the cuff with the velcro fastener. The cuff should fit snugly around the wrist. When relaxed, place your elbow on your chest with the wrist in the same height as your heart (see illustration 3).


Illustration 3

GB

50 Preparing to measure your blood pressure

 Do not attach the cuff over wounds as this can cause further injuries.

 Ensure that the cuff is not attached to an arm where arteries or veins are undergoing or have undergone medical treatment (e.g. shunt).

 For women who have had a breast amputated, do not attach the cuff to the wrist on the amputated side of the body.

 During the measurement, medical apparatus that are being used on the same arm at the same time, may malfunction.

 It can occur that there are differences between measurements at the wrist and the upper arm.

Therefore, it is recommended on occasions to take comparative measurements at the upper arm.

GB

Preparing to measure your blood pressure


51

 In cases of poor blood circulation or severe arteriosclerosis, it is not appropriate to take measurements at the wrist. In these instances the blood pressure should be taken at the upper arm.

 **When measuring, the wrist must be held at the same height as the heart!**

 The blood pressure measurement should not prevent the flow of blood for an unnecessary length of time (> 2 minutes). If the monitor fails to function correctly, remove the unit from the wrist.

 Excessively frequent measurements can cause injury by restricting blood flow.

52 Measuring your blood pressure

Your monitor has a memory with a capacity of 30 measurements.

➔ Start the measurement with the 'START' button

⚠ Now hold the arm absolutely still and do not talk.

All the items on the LCD display appear briefly to indicate that the monitor is now ready for use (see illustration 4).


Illustration 4

Measuring your blood pressure

53

The pump starts to inflate the cuff.

The display indicates the pressure of the cuff as it rises.

The unit inflates the cuff to a pressure of about 200 mmHg. The display indicates the pressure of the cuff as it rises. In case this pressure is not sufficient, the unit inflates the cuff by a further 50 mmHg.


Once inflated, the pump is switched off and the air in the cuff escapes automatically.

If not already the case, you should now keep your arm completely still and not speak.

As soon as the unit registers the first beats, the '♥' symbol blinks in time with the pulse rhythm.

Interrupting a measurement:

A measurement can be interrupted at any stage by pressing the 'START' button; the cuff deflates automatically.

54 Measurement display

After the end of the measurement, the inbuilt valve opens automatically and the cuff deflates rapidly.

The measured blood pressure values (systole, diastole, pulse) are displayed (see illustration 5).

Systole
mmHg

Diastole
mmHg

Pulse/min.

III. 5

The measurement is stored automatically (except for guest mode, see below). If the memory store is full (30 measurements), the **oldest** measurement is **deleted** and the current one stored.

Note on guest mode:

If a measurement value is not be kept, press „MEMORY“ while the measured values are still on display.

Measurement display

55

Invalid measurements with an 'Err' message (see page 58) are not stored.

If the symbol '❤️' appears after the measurement, it is recommended that the measurement is repeated, holding the arm completely still.

If the '❤️' symbol appears again even with the arm held still, this may be an indication of an irregular heartbeat. Please discuss this with your doctor at your next visit.

Blood pressure is a dynamic parameter and can be affected by the patient's position, e.g. sitting, standing, lying, moving, before or during the measurement as well as by his physical condition (stress, disease, etc.).


➔ If the measurement result is obviously wrong, repeat the measurement.

The unit switches itself off automatically after about 1 min. To repeat a measurement, press the 'START' button again.

Allow at least two minutes between measurements.

➔ If no further measurement is required, remove the unit from the wrist.

56 Completing the blood pressure record card

Retrieving the average value

➔ Enter each individual value in your record card.

stored measurements (see illustration 6).

➔ Recall the individual average value after 30 measurements by pressing the 'MEMORY' button for about 3 seconds while the unit is switched off. The average value of all stored values (systolic and diastolic) appears. The display <pulse> shows the number of

➔ Enter the average value in the field specifically reserved for this in your record card.

If there is no measurement in the memory, '00' is displayed for systole and diastole and 'A00' flashes for the number of measurements.

Illustration 6

Displaying individual stored values

57

- ➔ In order to bring back the last measurement stored, press the MEMORY button once more. The number of the measurements appears in the display. After three seconds the corresponding measurement value (systole, diastole and pulse) appears automatically.

Bild 7

- ➔ By repeated pressing of the MEMORY button, all stored measurements will be brought up in succession and appear in the display as described above.

If no button is pressed, the unit switches off automatically after 5 seconds.

- ➔ To start a new measurement, simply press „START“.
- ➔ **Please note:** if the batteries are removed the memory store is erased.


58 Error messages

If any problems occur during the measurement, the display will show an error message instead of a blood pressure reading.

Meaning of error messages:


Cause of Err message	Troubleshooting
 Inflation pressure in-adequate after three attempts to inflate.	Keep arm still
 No pulses were detected that could be evaluated.	Check the positioning of the cuff
 Error during inflation.	Cuff is possibly not tight enough

Error messages

"♥" Display	Troubleshooting
 Irregular puls or incorrect taking of blood pressure (i.e. movement during measurement)	Repeat measurement holding the arm absolutely still. If the "♥" symbol is displayed again, it can be that an irregular heart beat is present. Consult your GP.

60 Cleaning and Disinfection

Cleaning Monitor:

Use a soft, dry cloth to clean your monitor.

Cleaning Cuff:

Small stains on the cuff can be removed carefully with a proprietary washing-up liquid.


Disinfection Cuff:

For disinfectant wipes (at least 5 minutes exposure time) of the device and the cuff, we recommend the disinfectant antiseptic liquid (Schülke & Mayr). To disinfect the velcro of the cuff, we recommend spray disinfection. In particular, if the device is used on multiple users is to pay attention to regular cleaning and disinfections of the cuff.

Warranty

For details of the warranty conditions, please consult the warranty certificate supplied.

A claim under the warranty will only be accepted if the warranty certificate, completed and stamped by the dealer, is returned with the monitor.

Customer service

Warranty and repair work must be undertaken by trained and authorised personnel. Do not modify this equipment without authorization of the manufacturer.

Send the monitor, carefully packaged and with sufficient postage, to your authorised dealer or directly to:

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Disposal

61


Do not dispose of monitors and batteries in the domestic waste.

At the end of its lifespan, the monitor must be taken to a collection site for obsolete electronic items.

Please note that batteries and rechargeable batteries must be disposed of separately (local authority collection site).


Expected operational lifetime of the equipment: 10 years.


Expected operational lifetime of the cuff: 10.000 measurement cycles.

62 Technical data

Measurement principle:	oscillometric
Measurement range:	40 to 255 mmHg, 40 to 180 pulses/min.
Cuff pressure:	0 to 320 mmHg
 Memory store:	30 measurements
Display:	LCD
Operating conditions:	environmental temperature +10°C to +40°C relative humidity 15 to 85%
Storage conditions:	environmental temperature -10°C to +60°C relative humidity 15 to 85%


Technical data

63

Power supply:	DC 3 V (4 x 1.5 V Micro IEC LR 03 alkaline manganese batteries)
Typical battery life:	500 measurement cycles (depending on inflation pressure and frequency of use).
Battery check:	symbol display in LCD window
Weight:	150 g without batteries
Dimensions (WxHxD):	75 mm x 80 mm x 65 mm
Classification:	type BF ()
Protection against solid objects and water:	IP21

64 Technical data

**Clinical test
(DIN 58130):**

accuracy complies with the requirements
of EN 1060 Part 3

**Maximum deviation of cuff
pressure measurement:**

± 3 mmHg


**Maximum deviation of
pulse rate display:**

± 5 %

Applicable standards:

EN 1060, Part 1: 'Non-invasive sphygmomano-
meters – General requirements' and Part 3:
'Supplementary requirements for electro-
mechanical blood pressure measuring systems'

Calibration checks - Testing instructions

65

(for trained specialist staff only)

Calibration checks – every 2 years at the latest – must be performed by the following institutions or persons:

- Manufacturer
- Trade measurements authorities
- Persons who fulfil the requirements of the legislation governing the operation of medical devices.


A) Functional testing

Functional testing of the unit can only be performed on human subjects or using a suitable simulator.

B) Leak testing of the pressure circuit and error testing of the display

Please note:

- Prior to testing, separate the cuff from the unit (after having first removed the screws in the battery compartment). The air hose connecting plug (special plug) can now be inserted into the unit.

66 Calibration checks - Testing instructions

(for trained specialist staff only)


b) To perform leak testing of the pressure circuit, a tee connector must be used to connect the rolled-up cuff and a pump bulb to the pressure circuit.

c) If the pressure in measuring mode exceeds 320 mmHg, the quick-action air release valve triggers and opens the pressure circuit. If the pressure in testing mode exceeds 320 mmHg, the display blinks.

Testing

1. Remove batteries.
2. Hold START-button pressed down and insert batteries.
3. Release START-button.

The unit is now ready for testing. The displays of SYS, DIA and PUL now show the actual pressure.

Calibration checks - Testing instructions

67


(for trained specialist staff only)

4. Testing for divergence of Pressure indication and Density of Pressure circuit (setting time for cuff at least 60 seconds) can now be carried out in the usual manner.

C) Security

For security purposes the two parts of the housing are joined together with a security label.

5. Switch off unit and reconnect the cuff to the unit.

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BOSCH + SOHN GmbH u. Co. KG